



Figure 3: Picture of the setup. Dilatometer (left), Data logger (center) & VSP (right).

The five leads of the cell cable are plugged into dilatometer data logger box as shown in Fig. 4.





Figure 4: Cell cable connected to the data logger box.

The dilatometer and electrochemical data are recorded directly by EC-Lab® (no need to synchronize two software interfaces). The “External Devices” window has been configured as follows (Fig. 5):

- Select “Dilatometer” in the “Device Type” combo box.
- Select “ECD-nano-DL” in the “Device Name” drop box.

The settings of the device are already set by default.



Device Type: Dilatometer, Device Name: ECD-nano-DL

Analog OUT: Convert E/V to E/V, with 0 V = 0 V (max), 0 V = 0 V (min), E/V: 0.0

Analog IN 1: Convert E/V to Displacement/ μ m, with 10 V = 125 μ m (max), -10 V = -125 μ m (min)

Analog IN 2: Convert E/V to T/ $^{\circ}$ C, with 0.5 V = 100 $^{\circ}$ C (max), -0.125 V = -25 $^{\circ}$ C (min)


Figure 5: ECD-nano-DL Dilatometer configuration in the “External Devices” window.

Note:

ECD-nano-DL operated in 250 μ m (*i.e.* +/- 125 μ m) displacement range. Displacement range of 100 μ m is also available with the ECD-nano-DL.

III – MEASUREMENTS

The cycling of a LiCoO₂ cathode based battery is performed at 200 μ A in a potential window of 2.0 to 4.6 V. No floating period is set. The settings of the GCPL technique are shown in Fig. 6.



Set I to I_s = 0.200 mA vs. <None>

for at most t₁ = 1.00 h 0 mn 0.000 0 s

Limit Ewe > EM = 4.600 V

Record every dE₁ = 5.0 mV or dt₁ = 1.000 0 s

Hold EM for t_M = 0 h 0 mn 0.000 0 s

Limit ||| < I_M = 0.050 mA

Record every dQ = 1.000 mAh or dt_q = 1.000 0 s

Limit ΔQ > ΔQ_M = 0.000 mAh <=> Δx_M = 0.000

E Range = 0 V; 5 V Resolution = 100 μ V

I Range = 1 mA Bandwidth = 5

Figure 6: GCPL settings. Only the charge is shown in the screenshot.

Figure 7 shows the LiCoO_2 electrode dilation over seven successive cycles. As expected with the *in situ* XRD results [1], the dilation curve clearly shows an anomalous expansion during delithiation. Furthermore, the reported phase transitions close to the vertex points of the electrochemical cycle can be clearly seen. The phase transition at the upper vertex point at around 4.5 V fades out during cycling.




Figure 7: Electrode potential and dilation evolution during constant current cycling of a LiCoO_2 electrode. Bottom: focus on the phase transitions at the two vertex points and at the first discharge midpoint.

IV-CONCLUSION

This note presents how to connect the dilatometer and potentiostat/galvanostat together and perform an *in situ* dilation

measurement during battery cycling. The configuration of the EC-Lab[®] software is also de-scribed. It is noteworthy that similar measurements can be carried out on supercapacitors.

These coupled measurements are complementary to other characterizations such as X-ray measurements.

REFERENCES

- 1) Z. Chen, J. R. Dahn, *Electrochim. Acta*, 49(7) (2004) 1079.
- 2) J. N. Reimers, J. R. Dahn, *J. Electrochem. Soc.* 139(8) (1992) 2091.

Revised in 08/2019