

## Temperature calibration of a potentiostat board with an ITS

### I – INTRODUCTION


**Intermediate Temperature System (ITS)** is a compact temperature control unit based on the Peltier effect, used for the electrical characterization of materials. It is controlled by MT-Lab® software via an USB connection when used with the MTZ-35 impedance analyzer. The ITS is also controlled by a BioLogic potentiostat through its analog output using EC-Lab® software.

The ITS is temperature calibrated in the factory with the two provided PT1000 probes. For use with the MTZ-35 impedance analyzer no calibration is needed. However, for a use with a potentiostat, calibration may be needed because of a drift or an offset on the Analog out/Analog In of the user's potentiostat board.

This technical note details the procedure to follow for the calibration of the Analog out, the Analog In1 and the Analog In2 of a potentiostat board used for the temperature control and measurement of ITS.

### II– PRINCIPLE OF THE ITS CONTROL WITH A POTENTIOSTAT

The temperature of the ITS is controllable by a Potentiostat through the Analog out. Two temperature variables are measured: the ITS temperature termed here as " $T_{base}$ " and the cell temperature labelled as " $T_{cell}$ ". These two temperatures are acquired as two external voltages through the two Analog inputs of the potentiostat board. Figure 1 shows the " $T_{base}$ " and the " $T_{cell}$ " Temperature measured by the two PT-1000 provided with the ITS.



**Figure 1: Schematic showing the location of the PT-1000 probes and their corresponding temperatures.**

#### Note:


By default, two PT1000 are provided with the ITS. The first PT1000 is integrated into the ITS for measuring  $T_{base}$ . The second PT1000 is installed in the  $T_{cell}$  hole of CESH for measuring  $T_{cell}$ . An optional third PT1000 can be used with the ITS for measuring  $T_{sample}$ .

$T_{sample}$  can also be measured by EC-lab (Instead of  $T_{cell}$ ).

The calibration is carried out in two independent steps: a fast calibration of the Analog output and a second calibration of the analog inputs.

### III– CALIBRATION OF THE ANALOG OUTPUT

The default calibration parameters of the Analog out of the ITS are given in EC-Lab® software. They can be viewed by clicking on "External Devices" tab and selecting TCU in "Device Type" box and ITS in "Device Name" box (Fig. 2).



**Figure 2: Default settings of analog out for ITS.**

The calibration of the Analog out consists of plotting the graph “setpoint temperature” against a series of voltage values”. It is a linear relation between the temperature  $T$  and the voltage  $E$ :

$$T = a E + b \quad \text{where}$$

$T$  is the output signal (Temperature on ITS).  
 $a$  is the slope of the linear plot.  
 $E$  is the input signal (Voltage set in EC-Lab®).  
 $b$  is the offset temperature.

The calibration procedure of the Analog out is as follows:

- Connect the ITS to the potentiostat using the provided DB9-DB9 cable (Fig. 3).



Figure 3: Connection of the potentiostat board to an ITS using the DB9-DB9 cable.

- Open EC-Lab® software and connect it to the potentiostat.
- Load “OCV” technique (Fig. 4).
- Go to “External Devices” tab and select “Other” in “Device Type” and “Other” in “Device Name” box (Fig. 4).
- In Analog OUT section, select E/V variable in the “Convert” box and check it. Then set the values 10 V and 0 V as shown below.




Figure 4: Analog Out settings for calibration with ITS.

- Set a first setpoint voltage (for example 3.0 V) in the circled box on the Fig. 4 and press on the yellow key button to validate the setpoint voltage.
- Read and note the temperature setpoint displayed on the ITS screen.
- Set another setpoint voltage value and read the corresponding setpoint temperature displayed on the ITS screen.

The table I summarizes the voltage values and the corresponding setpoint temperatures obtained for the calibration of the Analog out.

Table I: Summary of the setpoint voltages set in the Analog OUT section and the corresponding setpoint temperatures displayed on the ITS screen.

E/V	ITS setpoint temperature/°C
3.0	- 9.9
4.0	23.4
5.0	56.8
6.0	90.1
7.0	123.4

The points ( $E, T$ ) are then used for plotting  $E$  vs.  $T$  linear graph (Fig. 5). The obtained linear relation is:  $E = 0.03 T + 3.297$  (Fig. 5).




Figure 5: Voltage vs. Temperature plot of Analog Out.

Once the linear relation has been determined, two points ( $T_1, E_1$ ) and ( $T_2, E_2$ ) have to be set in the EC-Lab® software. By default, the voltage at - 55 °C and 165 °C are usually used in EC-Lab® for the voltage-temperature conversion:

For  $T_1 = - 55$  °C;  $E_1 = 1.647$  V

For  $T_2 = 165\text{ }^\circ\text{C}$ ;  $E_2 = 8.247\text{ V}$

The values obtained here are the same as the default settings given in EC-Lab® (Fig. 2).

The setting with these two points has to be checked by using at least two temperatures.

**Note:**

The low and high temperature values set in Analog OUT section have to be lower than  $-40\text{ }^\circ\text{C}$  and higher than  $150\text{ }^\circ\text{C}$  respectively, as they are considered to be temperature limit values by EC-Lab® software.

A warm-up time of 30 mini is recommended before proceeding to the calibration. At this point, the calibration of Analog out can be performed quickly. It lasts less than 2 min.

**IV– CALIBRATION OF THE ANALOG INPUTS**

The  $T_{\text{base}}$  and the  $T_{\text{cell}}$  are acquired by Analog In1 and Analog In2 respectively. Analog In1 and Analog In2 operate in the voltage range  $\pm 10\text{ V}$ . The user needs at least two setpoint temperatures  $T_1$  and  $T_2$  to calibrate each Analog input.

**III-1 Step 1:**

The first step consists of creating two variables:  $E_1$  voltage for Analog In1 and  $E_2$  voltage for Analog In2. The user should proceed as follows:

- On EC-Lab® software, select “External Devices” tab.
- Click on the “Custom Units” button (Top right on Fig. 6). The window in Fig. 6 appears. Click on “Add” button.
- Type “ $E_1$ ” in the box “Name” and “V” in the box “Unit” (Fig. 6 left). Do the same thing for the variable “ $E_2$ ” (Fig. 6 right)




Figure 6: Addition of  $E_1$  and  $E_2$  voltage variables to the “External Devices” Tab.

**III-2 Step 2:**

Once the  $E_1$  and  $E_2$  variables have been created, the EC-Lab user can start the calibration of the Analog In1 and Analog In2 by proceeding as follows:

- The potentiostat has to be connected to ITS as shown in Fig. 3.
- Load Temperature control (TC) protocol on EC-Lab®.
- In “External Devices” tab, select TCU in “Device Type” box and ITS in “Device Name” box.
- Replace the  $T/^\circ\text{C}$  by  $E_1/\text{V}$  and  $E_2/\text{V}$  in the box in front of “Convert E/V” of Analog In1 & Analog In2 sections respectively and set the voltage values as shown in Fig. 7.




Figure 7: Analog In1 and Analog In2 settings for calibration.

- In Safety/Advanced Setting tab, check “smooth on” box with 10 points.
- In the Parameter settings tab, set the setpoint temperatures in the TC protocol (Fig. 8).

At least two temperature setpoints are needed to plot the linear Temp vs. Voltage graph for each Analog input. For a better accuracy three temperature setpoints were used here. Figure 8 shows the Temperature Control “TC” technique with three sequences corresponding to  $30\text{ }^\circ\text{C}$ ,  $60\text{ }^\circ\text{C}$  and  $90\text{ }^\circ\text{C}$  setpoint temperatures used for the calibration of the two Analog inputs.




Figure 8: TC protocol with three sequences used for the calibration of Analog In1 and Analog In2.

- Run the “TC” protocol and display the voltage versus time plots.

Once the steady state has been reached, the user has to read and note either the  $T_{base}$  and  $T_{cell}$  temperatures and their corresponding voltage values measured on EC-Lab®.

Figure 9 shows a picture of the ITS screen obtained with the setpoint temperatures 60.0 °C and 90.0 °C




Figure 9: Picture of the ITS screen showing the steady state values of  $T_{base}$  and  $T_{cell}$  obtained for the setpoint temperature of 60 °C (left) and 90 °C (right).

The voltage vs. time plots obtained with the three setpoint temperatures for analog In1 and Analog In2 are displayed on the Fig. 10.




Figure 10:  $E_1$  vs time and  $E_2$  vs. time plots obtained for a simultaneous calibration of Analog In1 & Analog In2.

The measured steady state temperatures ( $T_{base}$  and  $T_{cell}$ ) and voltages ( $E_1$  and  $E_2$ ) are summarized in the Table II below:


Table II: Summary of the  $T_{base}$  &  $T_{cell}$  temperatures and their corresponding voltages measured for each setpoint temperature.

ITS setpoint temperature/°C	$T_{base}$ /°C	$T_{cell}$ /°C	$E_1$ (Analog IN1) /mV	$E_2$ (Analog IN2) /mV
30.0	31.1	29.4	1754	1751
60.0	60.0	57.8	1945	1931
90.0	90.0	85.0	2136	2106

### Note

We strongly recommend you wait until the voltage  $E_1$  (and also  $E_2$ ) reaches a plateau. This means that  $T_{base}$  and  $T_{cell}$  have reached the equilibrium temperature. Indeed,  $T_{base}$  and  $T_{cell}$  displayed on ITS screen correspond to average temperature values whereas the  $T_{base}$  and  $T_{cell}$  acquired by EC-Lab® are instantaneous temperature values. When the steady state is reached, the instantaneous temperature values acquired by EC-Lab® are constant and are equal to the average temperature displayed on the ITS screen.

The Table II data were used for drawing the Voltage vs. Temperature plots of the calibration of Analog In1 and Analog In2. Figure 11 shows the obtained plots.






Figure 11: Voltage vs. Temperature plots due to the calibration of the Analog In1 (a) and Analog In2 (b).

The Voltage vs. Temperature plots were fitted to straight lines with the equations:  
 $E = 6.3768 T + 1562$  for Analog In1  
 $E = 6.3845 T + 1562$  for Analog In2

By default, the temperatures - 40 °C, 150 °C and the voltage at both temperatures are used as settings for Analog input in EC-Lab®. They are calculated using the calibration equations above.

Analog In1:

For  $T_1 = -40\text{ °C}$ ;  $E_1 = 1.307\text{ V}$   
 For  $T_2 = 150\text{ °C}$ ;  $E_2 = 2.519\text{ V}$

Analog In2:

$T_1 = -40\text{ °C}$ ;  $E_1 = 1.307\text{ V}$   
 $T_2 = 150\text{ °C}$ ;  $E_2 = 2.519\text{ V}$

These points are set in the Analog In1/In2 sections of the TCU device of EC-Lab®. The setting with these two points has to be checked for at least two temperatures.

## V– CONCLUSION

This note showed how a potentiostat user can calibrate analog out, analog In1 and Analog In2 of a potentiostat board used for the control of an ITS. This calibration improves the accuracy of temperature control and the accuracy of  $T_{\text{base}}$  and  $T_{\text{cell}}$  acquisitions and so leads to reliable temperature measurements.

## REFERENCES

- 1) Operating manual of ITS
- 2) [Technical note #2](#) CESH Sample holders.