

Applications of scanning probe electrochemistry– Sensors

SCAN-Lab

September 2021



Background




Background

9% of commercial scanning probe electrochemistry instruments are used in **sensors research.**

This document will further investigate the role of scanning probe electrochemistry in the field of sensors.

Publication Fields - All Techniques






Why is scanning probe electrochemistry applied in sensors?

Scanning probe electrochemistry is used in fundamental R&D studies of sensors to:

- Investigate the sensor analyte interaction
- Research novel sensor materials and structures
- Perform electrode modification
- Analyze ion selective membranes



What techniques are used?




Though there has been some application of **SKP** to sensor studies, **SECM** is by far the most popular scanning probe electrochemistry technique to be applied in this field.


Source: Analysis of scientific publications citing commercial instruments. Each research group was only counted once per technique.



What techniques to use for sensors?



SECM



SKP



Sensors



What are the research problems?

Bulk techniques are an average measure of the sensor system, and therefore do not provide insight on the **local properties of the sensor, e.g. **activity of the regions**.**

- Solution: Local measure of electrochemical activity

It is important to understand the effect and repeatability of **electrode modification.**

- Solution: Direct comparison of bare and modified electrodes

Traditional techniques to investigate biosensors require **sample pre-treatment, including **fluorescent labelling** limiting their applicability.**

- Solution: Measurement without pre-treatment

Understanding sensor systems requires knowledge of the **reactants, products, and ions involved in the system.**

- Solution: Chemical selectivity





Solution: Local measure of electrochemical activity

How this is met by scanning probe electrochemistry:

When sensor investigations are performed using bulk techniques the result is an average of the bulk sample. This can be difficult or impossible to interpret to determine the effect of local structure on activity towards an analyte.


In scanning probe electrochemistry **only the area under the probe is measured**, providing local data. The electrochemical characteristics measured by scanning probe electrochemistry relate to a sample's activity. By raster scanning the probe across the sample an x-y **map of activity** can be produced.



Solution: Direct comparison of bare and modified electrodes

How this is met by scanning probe electrochemistry:

Modification of sensor electrodes alters the electrochemical characteristics of the electrode. Using scanning probe electrochemistry the electrochemical activity of the sample can be **quantitatively mapped** before and after modification. These **maps can be directly compared** to see the local effects of modification on the sensor sample.






Solution: Measurement without pre-treatment

How this is met by scanning probe electrochemistry:


Many microscopies require the use of tags, or other pre-treatments which affect the system under study, and increase the required setup time. **Scanning probe electrochemistry uses characteristics inherent to the sample** to perform the measurement. SECM has the further advantage that it does not require electrical contact to the sample. All of this means that in scanning probe electrochemistry the **sample is measured as is**, without a lengthy preparation.



Solution: Chemical selectivity

How this is met by scanning probe electrochemistry:

Chemical selectivity is specifically met by **SECM**. In SECM the probe is **biased to interact selectively with a given electrochemically active species**. SECM can be used to perform quantitative measurements with chemical selectivity, with changes in current signal directly reflecting the concentration of the species of interest.





Conclusions



Summary

- Scanning probe electrochemistry is used to measure the **local electrochemical activity** of a sensor sample.
- Scanning probe electrochemistry offers a method of **directly comparing** the bare and modified electrodes.
- Scanning probe electrochemistry allows measurement **without pre-treatment**.
- Measurements with **chemical selectivity** are possible using scanning probe electrochemistry.
- In sensor research **SECM** is of the most interest, with **SKP** also used.



Why SCAN-Lab?

Using BioLogic's scanning probe electrochemistry instruments it is possible to locally measure sensor systems without the need for pre-treatment. BioLogic offers a wide range of techniques applicable to sensor research on a single, flexible, user-friendly platform. This expands the current and future research avenues available from a single instrument.



Learning Center Article

A series of Learning Center articles has been created to help determine the most appropriate technique for a given research problem. This includes an article dedicated to sensors research:

Scanning Probes & Sensor Research



Acronyms

- **SECM: Scanning ElectroChemical Microscopy**
- **SKP: Scanning Kelvin Probe**